

30µl Dispensing Tips


384-well disposable tips for high precision pipetting on BasePlate and Asset



30µl Dispensing Tips

Description

These 384-well format tips are designed to provide precise liquid dispensing of up to 30µl on the BasePlate and Asset automation platforms. They are supplied in a rack which enables the tips to be easily loaded into a standard metal carrier in the automated system.



Benefits

- Stepped body design minimizes volume displacement and ensures tip straightness for lower dead volumes
- Small internal tip orifice provides better small volume droplet formation, resulting in optimal accuracy and precision
- 100% inspection ensures every tip is free of defects
- Face seal design ensures consistent tip-to-tip height for precise low volume dispensing

Specifications

Tip diameter:	0.172"/4.37mm
Tip shoulder:	0.140"/3.56mm
Overall length:	1.384"/35.15mm
Inner orifice:	0.012"/0.305mm

Product Number

Cat No.	Description.
5351	30µl tips in traditional rack
5352	30µl tips in traditional rack - sterile

Minimum Order Quantity

1 case (10 racks of 384 tips)

How to Order

Please contact Thermo Scientific at:
US: +1-603-595-0505
Europe: +44-161-486-2110
or visit: www.thermo.com/matrix

UK and Mainland Europe
+44 (0) 1763 227 333

USA
+1 (302) 478 9060

Email
support@automationpartnership.com

