



Glass Solubilisation Tubes


Glass tubes for compounds and liquids offering
reliable identification and tracking



Glass Solubilisation Tubes

Description

Round bottom glass tubes for the solubilisation of compounds.



Benefits

- Uniquely bar coded for reliable identification and tracking
- Can be re-used
- Rack and sealing mats for secure storage and transfer are also available

Specification

- Manufactured from clear Borosilicate glass
- Round bottom tubes, 23.75mm diameter and 60mm length

Minimum Order Quantity

1 case (720 units)

How to Order

Contact Partnership Services quoting ordering code A-0020

Further Information

Racks and sealing mats are also available
See separate datasheet

UK and Mainland Europe
+44 (0) 1763 227 333

USA
+1 (302) 478 9060

Email
support@automationpartnership.com

