



# Microtube Caps

---



High quality, automation compatible microtube caps for proven reliable sealing



# Microtube Caps

## Description

High quality automation compatible caps for proven reliable sealing of microtubes. Their design enables easy handling by robotic equipment and a positive grip on the tube wall.



### Benefits

- Design ensures a positive fit on tube wall
- Compatible with tube handling automation
- Proven reliable seal with a range of microtubes

### Specification

- Injection moulded from medium density polyethylene (Esocrene)
- Compatible with standard 1.4ml microtubes
- Proven use for storing solubilized compounds in DMSO
- Compatible with The Automation Partnership's capper/de-capper automation

### Minimum Order Quantity

Available in cases of 100,000

### How to Order

Contact Partnership Services quoting order code 131-Z211-B

**UK and Mainland Europe**  
+44 (0) 1763 227 333

**USA**  
+1 (302) 478 9060

**Email**  
support@automationpartnership.com

