

# TAP Sealing Foil

---

Cost-effective sealing for a variety of plate types



# TAP Sealing Foil

## Description

Cost effective foil sealing material ideal for the automated heat sealers used in Compound Management and High Throughput Screening.



### Benefits

- Low pierce force required
- Seals both polypropylene and polyethylene plates and tubes
- Available in small quantities
- Resistant to DMSO when stored at below room temperature
- Compatible with The Automotive Partnership (TAP's) Picotube

### Specification


- Aluminium foil, 20 micron thickness
- 610 linear metres of foil per roll
- Compatible with most common sealers

### Minimum Order Quantity

1 reel of 610 linear metres

### How to Order

Contact Partnership Services quoting order code R-0061



**UK and Mainland Europe**

+44 (0) 1763 227 333

**USA**

+1 (302) 478 9060

**Email**

support@automationpartnership.com

