



Triple Flask Nunclon™ Δ

A high yield multi-layer cell culture flask compatible with
Select, Compact Select and Compact CellBase



Triple Flask

Description

Triple Flask has the same external dimensions as a T-175 flask but with three parallel growth surfaces providing a total cell culture area of 500cm². Existing cell types and processes can be readily transferred onto the Triple flask with the benefit of higher yield due to the increased surface area.



Benefits

- Increased cell yield, ideal for scale up
- Increased productivity being fully compatible within SelecT, CompacT SelecT and CompacT CellBase
- Increased storage capacity
- Reduced processing time

Specifications

- External footprint identical to a standard T-175 flask
- Robust. multi-layer construction
- Individually bar-coded for ease of tracking and auditing
- Nunclon™ Δ certified

Minimum Order Quantity

1 case

How to Order

Please contact Thermo Scientific at:
US: +1 (585) 899 7800
Europe: +45 4631 2000
or visit: www.nuncbrand.com

Product Number

Cat. No.	Growth Area (cm ²)	Surface Treatment	Neck Description	Cap	Qty/Pk	Qty/Cs
132920	500	Nunclon™ Δ	Straight	Filter	4	32

UK and Mainland Europe
+44 (0) 1763 227 333

USA
+1 (302) 478 9060

Email
support@automationpartnership.com

