

# Bulk Cell Output Tube Set

---

Tube set facilitating the rapid harvest and output of multiple flasks to an external suspension vessel

---

# Bulk Cell Output Tube Set

## Description

The tube set facilitates the rapid harvest and output of multiple flasks to an external suspension vessel, when used with the associated upgrade on Select and Compact.



## Benefits

Aseptic output of large quantities of cells held in suspension from Select and Compact

- Specially designed to collect the maximum cell culture output volumes of standard T175 flasks, Nunc Triple flask, Corning Expanded Surface Area flask and HYPERFlask
- Two collection vessels may be set up prior to each run allowing for the bulk output of two separate cell lines
- Minimises dead volumes in the flask

## Specification

- Superior grade silicone tubing manufactured under clean room (class 7 ISO 14644) conditions
- Gamma sterilised
- Two outputs may be fitted per Select
- Compatible with a variety of external suspension vessels and incubator configurations
- Designed to optimise the flow of cell culture into the collection vessel

## Minimum Order Quantity

1 box (10 packs) of individually sealed bags

Order code A-0024

## How to Order

Contact Partnership Services for further details

**UK and Mainland Europe**  
+44 (0) 1763 227333

**USA**  
+1 (302) 478 9060

**Email**  
support@automationpartnership.com

