



Large Volume Pipette Tip

50 μ l to 2500 μ l dispense volumes

Sterile 50 - 2500 μ l tip for use on the Cello liquid handler



Large Volume Pipette Tip

Description

A sterile 50-2500µl tip for use on the Cello liquid handler. Compatible with the 24-way cherry picking head and the 24-way fixed head.

Benefits

- Optimised tip design enables the same tip to be used for large volume liquid transfer into 6 and 24-well plates. Also provides precise control of small volume transfers into 96 and 384-well plates
- Narrow apex allows access into wells of 384-well plates
- Can dispense volumes as low as 50µl with 5% CV

Specification

- Transfer volume up to 2500µl per tip
- Made from polypropylene (compliant to United States Pharmacopeia XXII Class VI and European Pharmacopeia)
- Gamma sterilised (not suitable for autoclaving)
- Individually sealed in bag of 24 tips

Minimum Order Quantity

20 packs of 24 tips (480 tips)

How to Order

Contact Partnership Services quoting order code 032-L002-C



UK and Mainland Europe

+44 (0)1763 227 333

USA

+1 (302) 478 9060

Email

support@automationpartnership.com

